

Rabbit Anti-TREM1/AP Conjugated antibody

SL10306R-AP

Product Name	Anti-TREM1/AP
Chinese Name	碱性磷酸酶（AP）标记的髓系细胞触发受体 1 抗体
Alias	TREM 1; TREM-1; Triggering receptor expressed on monocytes 1; Triggering receptor expressed on myeloid cells 1; Triggering receptor TREM 1; Triggering receptor TREM1; Triggering-receptor TREM 1; Triggering-receptor TREM-1; Triggering-receptor TREM1; TREM1_HUMAN; CD354.
Research Area	Cell biology immunology Cell differentiation
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human(predicted:Rat) WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	24kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human TREM1
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail	background: This gene encodes a receptor belonging to the Ig superfamily that is expressed on myeloid cells. This protein amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as

increased surface expression of cell activation markers. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.[provided by RefSeq, Jun 2011].

Function:

Stimulates neutrophil and monocyte-mediated inflammatory responses. Triggers release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers. Amplifier of inflammatory responses that are triggered by bacterial and fungal infections and is a crucial mediator of septic shock.

Subunit:

Interacts with TYROBP/DAP12.

Subcellular Location:

Isoform 1: Cell membrane; Single-pass type I membrane protein (Potential).

Isoform 2: Secreted (Potential).

Tissue Specificity:

Highly expressed in adult liver, lung and spleen than in corresponding fetal tissue. Also expressed in the lymph node, placenta, spinal cord and heart tissues. Expression is more elevated in peripheral blood leukocytes than in the bone marrow and in normal cells than malignant cells. Expressed at low levels in the early development of the hematopoietic system and in the promonocytic stage and at high levels in mature monocytes. Strongly expressed in acute inflammatory lesions caused by bacteria and fungi. Isoform 2 was detected in the lung, liver and mature monocytes.

Post-translational modifications:

Glycosylated.

Similarity:

Contains 1 Ig-like V-type (immunoglobulin-like) domain.

Database links:

[Entrez Gene: 54210](#) Human

[Entrez Gene: 58217](#) Mouse

[Omim: 605085](#) Human

[SwissProt: Q9NP99](#) Human

[SwissProt: Q9JKE2](#) Mouse



[Unigene: 283022](#) Human

[Unigene: 248352](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.