

Rabbit Anti-2-Naphthoic acid antibody

SL10418R

Product Name 2-Naphthoic acid

Chinese Name 2-萘甲酸抗体

Alias 2-Naphthalenecarboxylic acid; β -Naphthoic acid; Isonaphthoic acid; 2-Carboxynaphthalene;
2-Maythic acid; 2-Naphthoic acid; Naphthalene- β -carboxylic acid; NSC 59901;
2-Naphthylcarboxylic acid

Product Type Small molecule anti

Immunogen Species Rabbit

Clonality Polyclonal

React Species Species independent

Applications IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated 2-Naphthoic acid

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution Species independent1M TBS(pH7.4) with 1% BSA, Species independent3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

2-Naphthoic acid (NPA) is a noncompetitive N-methyl-D-aspartate (NMDA) receptor inhibitor. The fluorescence spectra and electronic absorption of 2-naphthoic acid was studied.

Product Detail

SWISS:
N/A



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

CAS:
93-09-4