

## Rabbit Anti-parathyroid hormone antibody

SL10488R

**Product Name** parathyroid hormone

**Chinese Name** 甲状旁腺激素抗体

**Alias** hPTH; Parathormone; Parathyrin; Parathyroid hormone 1; Parathyroid hormone; PTH; PTH1; PTH1 receptor; PTH1R; PTHR; PTHR1; PTHY\_HUMAN.

**Immunogen Species** Rabbit

**Clonality** Polyclonal

**React Species** (predicted: Human, Mouse, Rat, )

**Applications** IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000  
(Paraffin sections need antigen repair)  
not yet tested in other applications.  
optimal dilutions/concentrations should be determined by the end user.

**Theoretical molecular weight** 13kDa

**Cellular localization** Secretory protein

**Form** Liquid

**Concentration** 1mg/ml

**immunogen** KLH conjugated synthetic peptide derived from human parathyroid hormone: 11-80/115

**Lsotype** IgG

**Purification** affinity purified by Protein A

**Buffer Solution** 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

**Storage** Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

**PubMed** [PubMed](#)

**Product Detail**