

## Rabbit Anti-BACE2/AP Conjugated antibody

SL10494R-AP

<b>Product Name</b>	Anti-BACE2/AP
<b>Chinese Name</b>	碱性磷酸酶 (AP) 标记的 $\beta$ 淀粉样前体蛋白裂解酶 2 抗体
<b>Alias</b>	ALP 56; ALP-56; ALP56; ALP56; Asp 1; ASP1; ASP1; ASP21; Aspartic Like Protease 56kDa; Aspartic Like Protease 56kDa; Aspartic-like protease 56 kDa; Aspartyl protease 1; Aspartyl protease 1; BACE 2; BACE 2; Bace2; BACE2_HUMAN; Beta Secretase 2; Beta Secretase 2; Beta Site Amyloid Beta A4 Precursor Protein Cleaving Enzyme 2; Beta Site Amyloid Beta A4 Precursor Protein Cleaving Enzyme 2; Beta-secretase 2; Beta-site amyloid precursor protein cleaving enzyme 2; Beta-site APP cleaving enzyme 2; Beta-site APP cleaving enzyme 2; Down region aspartic protease; Down region aspartic protease; Downs Syndrome Region Aspartic Protease; Downs Syndrome Region Aspartic Protease; DRAP; Memapsin 1; Memapsin1; Memapsin-1; Membrane associated aspartic protease 1; Membrane-associated aspartic protease 1; Theta-secretase.
<b>Research Area</b>	
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Human,Mouse,Rat,Dog,Pig,Cow,Sheep) IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:500-5000
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	50kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human BACE2
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH

7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

**background:**

Cerebral deposition of amyloid beta peptide is an early and critical feature of Alzheimer's disease and a frequent complication of Down syndrome. Amyloid beta peptide is generated by proteolytic cleavage of amyloid precursor protein by 2 proteases, one of which is the protein encoded by this gene. This gene localizes to the 'Down critical region' of chromosome 21. The encoded protein, a member of the peptidase A1 protein family, is a type I integral membrane glycoprotein and aspartic protease. Three transcript variants encoding different isoforms have been described for this gene. [provided by RefSeq, Jul 2008].

**Function:**

Responsible for the proteolytic processing of the amyloid precursor protein (APP). Cleaves APP, between residues 690 and 691, leading to the generation and extracellular release of beta-cleaved soluble APP, and a corresponding cell-associated C-terminal fragment which is later released by gamma-secretase. It has also been shown that it can cleave APP between residues 671 and 672.

**Subunit:**

**Product Detail**

Monomer. Interacts with RTN3 and RTN4.

**Subcellular Location:**

Membrane. Golgi apparatus. Endoplasmic reticulum. Endosome. Cell surface.

**Tissue Specificity:**

Brain. Present in neurons within the hippocampus, frontal cortex and temporal cortex (at protein level). Expressed at low levels in most peripheral tissues and at higher levels in colon, kidney, pancreas, placenta, prostate, stomach and trachea. Expressed at low levels in the brain. Found in spinal cord, medulla oblongata, substantia nigra and locus coeruleus. Expressed in the ductal epithelium of both normal and malignant prostate.

**Post-translational modifications:**

Undergoes autoproteolytic cleavage.  
Glycosylated.

**Similarity:**

Belongs to the peptidase A1 family.

**Database links:**



[Entrez Gene: 25825](#) Human

[Entrez Gene: 56175](#) Mouse

[Entrez Gene: 288227](#) Rat

[Omim: 605668](#) Human

[SwissProt: Q9Y5Z0](#) Human

[SwissProt: Q9JL18](#) Mouse

[SwissProt: Q6IE75](#) Rat

[Unigene: 529408](#) Human

[Unigene: 97885](#) Mouse

[Unigene: 38581](#) Rat

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.