



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

## Rabbit Anti-Listeriolysin antibody

SL10499R

<b>Product Name</b>	Listeriolysin
<b>Chinese Name</b>	李斯特菌溶胞素抗体
<b>Alias</b>	Listeriolysin; hly; hlyA; lisA; LLO; Thiol activated cytolysin; listeriolysin precursor.
<b>Research Area</b>	Microbiology Bacteria and viruses
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Listeria monocytogenes) IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000 (Paraffin sections need antigen repair)
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Theoretical molecular weight</b>	55kDa
<b>Cellular localization</b>	The cell membrane Secretory protein
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from Listeria monocytogenes Listeriolysin: 251-350/529
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	(predicted:Listeria monocytogenes)1M TBS(pH7.4) with 1% BSA, (predicted:Listeria monocytogenes)3% Proclin300 and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
<b>PubMed Product</b>	<a href="#">PubMed</a>



SunLong Biotech Co.,LTD  
Tel: 0086-571-56623320 Fax:0086-571-56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

---

## **Detail**