

Rabbit Anti-DIM-7 antibody

SL10513R

Product Name DIM-7

Chinese Name DIM-7 蛋白抗体

Alias moleskin; D-Importin 7/RanBP7; DIM-7; Dmel_CG7935; 61; CG7935; CIP-61; Dim-7; Dim7; DmRanBP7; Dmel\CG7935; Imp7; Imp7/8; Msk; S(ls)2; dim-7; dim7.

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted: Fruit Fly,)

Applications IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 127kDa

Cellular localization The nucleus

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from fruit fly DIM-7: 401-500/1049

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail