

Rabbit Anti-blue-sensitive visual pigment antibody

SL10534R

Product Name blue-sensitive visual pigment

Chinese Name 蓝色敏感的视色素 (厚皮螺旋体) 抗体

Alias blue-sensitive visual pigment [Heliconius pachinus]; 蓝敏视觉色素; 厚皮螺旋体;

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted:Heliconius pachinus)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA
(Paraffin sections need antigen repair)

Applications not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 42kDa

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from Heliconius pachinus blue-sensitive visual pigment

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution (predicted:Heliconius pachinus)1M TBS(pH7.4) with 1% BSA, (predicted:Heliconius pachinus) and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail