



Rabbit Anti-Long wavelength-sensitive opsin antibody

SL10535R

Product Name Long wavelength-sensitive opsin

Chinese Name 长波长敏感的视蛋白抗体

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted:Heliconius pachinus)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA

Applications (Paraffin sections need antigen repair)

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 42kDa

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from Heliconius pachinusLong wavelength-sensitive 31-130/381

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution (predicted:Heliconius pachinus)1M TBS(pH7.4) with 1% BSA, (predicted:Heliconius pachinus) and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail