



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

Rabbit Anti-Interferon alpha 8 antibody

SL10573R

Product Name Interferon alpha 8

Chinese Name Interferon α 8 抗体

Alias IFN Alpha 8; IFNA8; IFN-alpha-8; IFN-alpha 8; IFNA8; IFN-alphaB; Interferon alpha B; Interferon alpha B2; interferon alpha type 201; interferon alpha-B'; interferon, alpha 8; LeIF B; OTTHUMP00000021146; IFNA8_HUMAN.

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted: Human,)

Applications IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)
not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 18kDa

Cellular localization Secretory protein

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from human Interferon alpha 8

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)

Product Detail



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com
