

## Rabbit Anti-Ear5 antibody

SL10605R

**Product Name** Ear5

**Chinese Name** 嗜酸性粒细胞阳离子型核糖核酸酶 5 抗体

**Alias** eosinophil cationic-type ribonuclease 5 precursor; ECP5\_MOUSE; Eosinophil cationic-type ribonuclease 5; MR-5; Ear5; Ear 5; Ear-5; mR5.

**Research Area** Cell biology immunology

**Immunogen Species** Rabbit

**Clonality** Polyclonal

**React Species** Mouse(predicted:Human,Rat)

WB=1:500-2000,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,Flow-Cyt=1ug/Test  
(Paraffin sections need antigen repair)

**Applications** not yet tested in other applications.  
optimal dilutions/concentrations should be determined by the end user.

**Theoretical molecular weight** 14kDa

**Form** Liquid

**Concentration** 1mg/ml

**immunogen** KLH conjugated synthetic peptide derived from mouse Ear5: 101-155/155

**Lsotype** IgG

**Purification** affinity purified by Protein A

**Buffer Solution** 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

**Storage** Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

**PubMed** [PubMed](#)

**Product Detail**