

## Rabbit Anti-MyT1L antibody

SL10609R

**Product Name** MyT1L

**Chinese Name** 神经锌指转录因子 1 抗体

**Alias** Myelin transcription factor 1-like protein; MyT1-L; MyT1L; Myt1l; MYT1L\_HUMAN; neural zinc finger transcription factor 1.

**Research Area** Cell biology immunology Neurobiology

**Immunogen Species** Rabbit

**Clonality** Polyclonal

**React Species** (predicted: Human, Mouse, Rat, Chicken, Dog, Pig, Cow, Horse, Sheep, )  
IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000

**Applications** (Paraffin sections need antigen repair )  
not yet tested in other applications.  
optimal dilutions/concentrations should be determined by the end user.

**Theoretical molecular weight** 133kDa

**Cellular localization** The nucleus

**Form** Liquid

**Concentration** 1mg/ml

**immunogen** KLH conjugated synthetic peptide derived from human MyT1L: 901-1000/1186

**Lsotype** IgG

**Purification** affinity purified by Protein A

**Buffer Solution** 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

**Storage** Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

**PubMed** [PubMed](#)

**Product**



SunLong Biotech Co.,LTD  
Tel: 0086-571-56623320 Fax:0086-571-56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

---

## **Detail**