

Rabbit Anti-ZmCCD8 antibody

SL10630R

Product Name	ZmCCD8
Chinese Name	类胡萝卜素裂解双加氧酶 8 抗体
Alias	Carotenoid cleavage dioxygenase 8;
Product Type	Botany 抗体
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Zea mays) IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500,ELISA=1:5000-10000 (Paraffin sections need antigen repair)
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	63kDa
Form	Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human ZmCCD8: 441-540/572
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:Zea mays)1M TBS(pH7.4) with 1% BSA, (predicted:Zea mays)3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	