

Rabbit Anti-Ultraviolet wavelength-sensitive opsin/Cy5.5 Conjugated antibody

SL10729R-Cy5. 5

Product Name	Anti-Ultraviolet wavelength-sensitive opsin/Cy5.5
Chinese Name	Cy5.5 标记的紫外波长敏感的视蛋白抗体
Alias	UVRh.
Research Area	
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Helicoverpa armigera) ICC/IF=1:50-200,IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	42kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human Ultraviolet wavelength-sensitive opsin
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Important Note: This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.