

Rabbit Anti-TMEM87A/AP Conjugated antibody

SL12666R-AP

Product Name	Anti-TMEM87A/AP
Chinese Name	碱性磷酸酶（AP）标记的 Transmembrane protein87A 抗体
Alias	A930025J12Rik; PSEC0094; RGD1309417; RP23-348L16.13-003; TM87A_HUMAN; tmem87a; Transmembrane protein 87A.
Research Area	Cell biology immunology Transmembrane protein
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Mouse(predicted:Human,Rat,Chicken,Dog,Pig,Cow,Horse,Rabbit,Sheep) WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	61kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human TMEM87A
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
	Subcellular Location: Membrane.
Product Detail	Similarity: Belongs to the LU7TM family. TMEM87 subfamily.

Database links:

[Entrez Gene: 25963](#) Human

[SwissProt: Q8NBN3](#) Human

[Unigene: 511138](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.