

Rabbit Anti-SAM68/Cy5 Conjugated antibody

SL13699R-Cy5

Product Name	Anti-SAM68/Cy5
Chinese Name	Cy5 标记的 SRC 有丝分裂相关蛋白 68 抗体
Alias	FLJ34027; GAP associated tyrosine phosphoprotein p62; GAP-associated tyrosine phosphoprotein p62; KH domain containing RNA binding signal transduction associated 1; KH domain-containing; KHDR1_HUMAN; KHDRBS1; p21 Ras GTPase activating protein associated p62; p21 Ras GTPase-activating protein-associated p62; p62; p68; RNA-binding; Sam68; signal transduction-associated protein 1; Src associated in mitosis 68 kDa protein; Src-associated in mitosis 68 kDa protein.
Research Area	Cell biology Signal transduction transcriptional regulatory factor Epigenetics
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Rat(predicted:Human,Mouse,Chicken,Dog,Pig,Cow,Horse,Rabbit,Sheep) IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	68kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human SAM68
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 癆 for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20 癆. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 癆.
Product Detail	background: Sam 68 is a protein that is phosphorylated on tyrosine and functions as a

substrate for Src family tyrosine kinases during mitosis. Sam 68 also associates with several SH2 and SH3 domain-containing signaling proteins, such as GRB2 and PLC γ . Originally cloned as Ras GAP-associated p62, further investigations have shown that Sam 68 and Ras GAP-associated p62 are not antigenically related, nor are they encoded by the same gene. Like Sam 68, the Sam 68-like mammalian proteins, SLM-1 and SLM-2, demonstrate RNA binding activity. Also like Sam 68, SLM-1 is tyrosine phosphorylated and functions as an adapter protein for signaling molecules, including GRB2, PLC γ , Fyn and Ras GAP. SLM-2 is not tyrosine phosphorylated, nor does it appear to associate with GRB2, PLC γ , Fyn or Ras GAP, indicating that SLM-2 may not be an adapter protein for these proteins.

Function:

Recruited and tyrosine phosphorylated by several receptor systems, for example the T-cell, leptin and insulin receptors. Once phosphorylated, functions as an adapter protein in signal transduction cascades by binding to SH2 and SH3 domain-containing proteins. Role in G2-M progression in the cell cycle. Represses CBP-dependent transcriptional activation apparently by competing with other nuclear factors for binding to CBP. Also acts as a putative regulator of mRNA stability and/or translation rates and mediates mRNA nuclear export. Isoform 3, which is expressed in growth-arrested cells only, inhibits S phase.

Subcellular Location:

Nucleus. Membrane.

Tissue Specificity:

Ubiquitously expressed in all tissue examined. Isoform 1 is expressed at lower levels in brain, skeletal muscle, and liver whereas isoform 3 is intensified in skeletal muscle and in liver.

Post-translational modifications:

Tyrosine phosphorylated by several non-receptor tyrosine kinases, for example LCK, FYN and JAK3. Negatively correlates with ability to bind RNA but required for many interactions with proteins.

Acetylated. Positively correlates with ability to bind RNA.

Arginine methylation is required for nuclear localization. Also can affect interaction with other proteins. Inhibits interaction with Src-like SH3 domains, but not interaction with WW domains of WBP4/FBP21 AND FNBP4/FBP30. Arg-291, Arg-331 and Arg-346 are found to be also dimethylated, probably to asymmetric dimethylarginine.

Similarity:

Belongs to the KHDRBS family.

Contains 1 KH domain.

Database links:

[Entrez Gene: 374071](#) Chicken

[Entrez Gene: 10657](#) Human

[Entrez Gene: 20218](#) Mouse

[Entrez Gene: 117268](#) Rat

[Omim: 602489](#) Human

[SwissProt: Q8UUW7](#) Chicken

[SwissProt: Q07666](#) Human

[SwissProt: Q3U8T3](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.