

Rabbit Anti-DDX60L antibody

SL14236R

Product Name DDX60L

Chinese Name ATP 依赖 RNA 解旋酶 DDX60 样蛋白抗体

Alias DDX60L; DDX6L_HUMAN; DEAD (Asp Glu Ala Asp) box polypeptide 60 like; DEAD box protein 60-like; DKFZp781D1175; FLJ13468; FLJ31033; FLJ39050; Probable ATP-dependent RNA helicase DDX60-like.

Research Area Cell biology The new supersedes the old

Immunogen Species Rabbit

Clonality Polyclonal

React Species (predicted: Human,)

IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000
(Paraffin sections need antigen repair)

Applications not yet tested in other applications.
optimal dilutions/concentrations should be determined by the end user.

Theoretical molecular weight 198kDa

Cellular localization cytoplasmic

Form Liquid

Concentration 1mg/ml

immunogen KLH conjugated synthetic peptide derived from human DDX60L: 501-600/1706

Lsotype IgG

Purification affinity purified by Protein A

Buffer Solution 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

Storage Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

Attention This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

PubMed [PubMed](#)



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

**Product
Detail**