

## Rabbit Anti-DHRS4L1 antibody

SL14309R

**Product Name** DHRS4L1

**Chinese Name** 短链脱氢酶/还原酶家族 4 样 1 蛋白抗体

**Alias** DHRS4L1; DR4L1\_HUMAN; Putative dehydrogenase/reductase SDR family member 4-like 2.

**Research Area** Cell biology immunology

**Immunogen Species** Rabbit

**Clonality** Polyclonal

**React Species** (predicted: Human, Mouse, Rat, Dog, Pig, Horse, Rabbit, Sheep, )  
IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500,ELISA=1:5000-10000  
(Paraffin sections need antigen repair)  
**Applications** not yet tested in other applications.  
optimal dilutions/concentrations should be determined by the end user.

**Theoretical molecular weight** 31kDa

**Cellular localization** cytoplasmic

**Form** Liquid

**Concentration** 1mg/ml

**immunogen** KLH conjugated synthetic peptide derived from human DHRS4L1: 151-250/281

**Lsotype** IgG

**Purification** affinity purified by Protein A

**Buffer Solution** 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

**Storage** Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

**PubMed** [PubMed](#)

**Product**



SunLong Biotech Co.,LTD  
Tel: 0086-571-56623320 Fax:0086-571-56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

---

## **Detail**