

Rabbit Anti-CANP/AP Conjugated antibody

SL14694R-AP

Product Name	Anti-CANP/AP
Chinese Name	碱性磷酸酶 (AP) 标记的癌相关核蛋白抗体
Alias	FAM111B; F111B_HUMAN; Cancer associated nucleoprotein; CANP; Protein FAM111B.
Research Area	Tumour Cell biology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Human) IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:500-5000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	85kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human CANP
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	background: FAM111B belongs to the FAM111 family. There are two named isoforms.
Product Detail	Similarity: Belongs to the FAM111 family. Database links:



[Entrez Gene: 374393](#) Human

[SwissProt: Q6SJ93](#) Human

[Unigene: 186579](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.