

Rabbit Anti-FAM47A/PE Conjugated antibody

SL14996R-PE

Product Name	Anti-FAM47A/PE
Chinese Name	PE 标记的 FAM47A 蛋白抗体
Alias	Family with sequence similarity 47 member A; FA47A_HUMAN; Hypothetical protein LOC158724; MGC27003.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	90kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human FAM47A
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	background: FAM47A belongs to the FAM47 family. The exact function of FAM47A remains unknown.
Product Detail	Database links: Entrez Gene: 158724 Human



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

[SwissProt: Q5JRC9](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.