

Rabbit Anti-FAM47E/AP Conjugated antibody

SL14998R-AP

Product Name	Anti-FAM47E/AP
Chinese Name	碱性磷酸酶 (AP) 标记的 FAM47E 蛋白抗体
Alias	FA47E_HUMAN; FAM47E; family with sequence similarity 47, member E; Protein FAM47E.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Human) IHC-P=1:100-500,IHC-F=1:100-500,ELISA=1:500-5000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	46kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human FAM47E
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
	Database links:
	Entrez Gene: 100129583 Human
Product Detail	SwissProt: Q6ZV65 Human
	Unigene: 720509 Human



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.