

## Rabbit Anti-C3Orf34/Cy5 Conjugated antibody

SL15175R-Cy5

<b>Product Name</b>	Anti-C3Orf34/Cy5
<b>Chinese Name</b>	Cy5 标记的 3 号染色体开放阅读框 34 抗体
<b>Alias</b>	Chromosome 3 open reading frame 34; HSD5; LOC84984; MGC14126; CEP19_HUMAN.
<b>Research Area</b>	Tumour Cell biology immunology
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	Mouse(predicted:Human,Rat,Chicken,Dog,Pig,Cow,Horse) IF=1:100-500
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	19kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration immunogen</b>	1mg/ml KLH conjugated synthetic peptide derived from human C3Orf34
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Product Detail</b>	<b>background:</b> C3orf34 (chromosome 3 open reading frame 34), also known as MGC14126, is a 163 amino acid protein encoded by a gene that maps to human chromosome 3q29. Chromosome 3 is made up of approximately 214 million bases encoding over 1,100 genes. Notably, there is a chemokine receptor gene cluster and a variety of human cancer related loci on chromosome 3. Particular

regions of the chromosome 3 short arm are deleted in many types of cancer cells. Key tumor suppressing genes on chromosome 3 encode apoptosis mediator RASSF1, cell migration regulator HYAL1 and angiogenesis suppressor SEMA3B. Marfan Syndrome, porphyria, von Hippel-Lindau syndrome, osteogenesis imperfecta and Charcot-Marie-Tooth disease are a few of the numerous genetic diseases associated with chromosome 3.

**Subcellular Location:**

Cytoplasm, cytoskeleton, centrosome, centriole. Cytoplasm, cytoskeleton, spindle pole. Note=Associates with the mother centriole in early interphase. Localizes to spindle poles during mitosis, and to distinct foci oriented.

**Similarity:**

Belongs to the CEP19 family.

**Database links:**

[Entrez Gene: 84984](#) Human

[SwissProt: Q96LK0](#) Human

[SwissProt: Q9CQA8](#) Mouse

[SwissProt: Q9QZX9](#) Rat

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.