

Rabbit Anti-C6orf223/AP Conjugated antibody

SL15242R-AP

Product Name	Anti-C6orf223/AP
Chinese Name	碱性磷酸酶（AP）标记的 6 号染色体开放阅读框 223 抗体
Alias	Chromosome 6 open reading frame 223; Hypothetical protein LOC221416; MGC45491; Uncharacterized protein C6orf223; CF223_HUMAN.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Human) ELISA=1:500-5000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	26kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human c6orf223
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
	Database links:
	Entrez Gene: 221416 Human
Product Detail	SwissProt: Q8N319 Human



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.