

Rabbit Anti-HEATR5B/PE Conjugated antibody

SL15441R-PE

Product Name	Anti-HEATR5B/PE
Chinese Name	PE 标记的 HEAT 重复内含蛋白 5B 蛋白抗体
Alias	HEAT repeat-containing protein 5B; heatr5b; HTR5B_HUMAN; KIAA1414.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Human,Mouse,Rat,Chicken,Dog,Pig,Cow,Horse,Rabbit,Sheep) ICC/IF=1:50-200,IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	224kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human HEATR5B
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Similarity: Belongs to the HEATR5 family. Contains 3 HEAT repeats.

Database links:



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

[Entrez Gene: 54497](#) Human

[SwissProt: Q9P2D3](#) Human

[Unigene: 468186](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.