



Rabbit Anti-HEATR8/AP Conjugated antibody

SL15443R-AP

| | |
|--------------------------|--|
| Product Name | Anti-HEATR8/AP |
| Chinese Name | 碱性磷酸酶（AP）标记的 HEAT 重复内含蛋白 8 蛋白抗体 |
| Alias | C1orf175; HEAT repeat-containing protein 8; MROH7_HUMAN; Heatr8; |
| Research Area | Cell biology immunology |
| Immunogen Species | Rabbit |
| Clonality | Polyclonal |
| React Species | (predicted:Human,Mouse,Rat,Dog,Pig,Rabbit) ELISA=1:500-5000 |
| Applications | not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Molecular weight | 146kDa |
| Form | Lyophilized or Liquid |
| Concentration | 1mg/ml |
| immunogen | KLH conjugated synthetic peptide derived from human HEATR8 |
| Lsotype | IgG |
| Purification | affinity purified by Protein A |
| Storage Buffer | 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C. |
| Storage | Subcellular Location: Membrane; Multi-pass membrane protein (Potential). |
| Product Detail | Similarity: Contains 4 HEAT repeats. Database links: |



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

[Entrez Gene: 374977](#) Human

[SwissProt: Q68CQ1](#) Human

[Unigene: 412482](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.