

Rabbit Anti-KCTD19/AP Conjugated antibody

SL16930R-AP

Product Name	Anti-KCTD19/AP
Chinese Name	碱性磷酸酶（AP）标记的钾离子通道多聚体结构域蛋白 19 抗体
Alias	BTB/POZ domain-containing protein KCTD19; FLJ40162; KCD19_HUMAN; KCTD19; potassium channel tetramerisation domain containing 19.
Research Area	Cell biology immunology Neurobiology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Mouse,Rat(predicted:Human,Dog,Pig,Cow,Horse,Rabbit,Sheep)
Applications	IHC-P=1:100-500,IHC-F=1:100-500 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	105kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human KCTD19
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	Similarity: Contains 2 BTB (POZ) domains.
Product Detail	Database links: Entrez Gene: 146212 Human



[Entrez Gene: 279499](#) Mouse

[Entrez Gene: 291965](#) Rat

[SwissProt: Q17RG1](#) Human

[SwissProt: Q562E2](#) Mouse

[Unigene: 299127](#) Human

[Unigene: 67628](#) Mouse

[Unigene: 144770](#) Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.