

Rabbit Anti-SPANXN2/Cy5 Conjugated antibody

SL17620R-Cy5

Product Name	Anti-SPANXN2/Cy5
Chinese Name	Cy5 标记的核相关精子蛋白 N2 抗体
Alias	Nuclear-associated protein SPAN-Xn2; SPANX family member N2; SPANX-N2; SPANXN2; Sperm protein associated with the nucleus on the X chromosome N2; SPXN2_HUMAN.
Research Area	Cell biology Developmental biology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	(predicted:Human,Chimpanzee) ICC/IF=1:50-200,IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	20kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human SPANXN2
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Similarity: Belongs to the SPAN-X family. Database links:



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

[Entrez Gene: 494119](#) Human

[SwissProt: Q5MJ10](#) Human

[Unigene: 535083](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.