

Rabbit Anti-IFNA17/Cy5.5 Conjugated antibody

SL18123R-Cy5. 5

Product Name	Anti-IFNA17/Cy5.5
Chinese Name	Cy5.5 标记的 Interferon α 17 抗体
Alias	IFN-alpha-17; IFN17_HUMAN; IFNA17; INF-alpha-17; interferon alpha 17 precursor; interferon alpha 88; interferon alpha I'; interferon alpha T; interferon alpha WA; Interferon alpha-17; Interferon alpha-88; Interferon alpha-I''; Interferon alpha-T; interferon, alpha 17; LeIF 1; LeIF I.
Research Area	Tumour Cell biology immunology Signal transduction
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human
Applications	IF=1:100-500 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	19kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human IFNA17
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail	background: Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.



Subcellular Location:

Secreted.

Similarity:

Belongs to the alpha/beta interferon family.

Database links:

[Entrez Gene: 3451](#) Human

[Omim: 147583](#) Human

[SwissProt: P01571](#) Human

[Unigene: 282276](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.