

Rabbit Anti-TNRC18/AP Conjugated antibody

SL18159R-AP

Product Name	Anti-TNRC18/AP
Chinese Name	碱性磷酸酶（AP）标记的 TNRC18 蛋白抗体
Alias	Long CAG trinucleotide repeat-containing gene 79 protein; TNC18_HUMAN; TNRC18; Trinucleotide repeat containing 18; Trinucleotide repeat-containing gene 18 protein.
Research Area	Cell biology transcriptional regulatory factor Binding protein Epigenetics
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Mouse,Rat(predicted:Human) IHC-P=1:100-500,IHC-F=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	315kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human TNRC18
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	Similarity: Contains 1 BAH domain.
Product Detail	Database links: Entrez Gene: 84629 Human



[Entrez Gene: 231861](#) Mouse

[SwissProt: O15417](#) Human

[SwissProt: Q80WC3](#) Mouse

[Unigene: 721291](#) Human

[Unigene: 354649](#) Mouse

[Unigene: 480618](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.