

Rabbit Anti-ZDHHC6/Cy5.5 Conjugated antibody

SL18485R-Cy5. 5

Product Name	Anti-ZDHHC6/Cy5.5
Chinese Name	Cy5.5 标记的 ZDHHC6 蛋白抗体
Alias	DHHC-6; FLJ21952; Probable palmitoyltransferase ZDHHC6; Transmembrane protein H4; ZDHC6_HUMAN; Zdhhc6; Zinc finger DHHC domain-containing protein 6; Zinc finger protein 376; ZNF376.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human,Mouse,Rat(predicted:Pig,Horse) ICC=1:50-200 IF=1:50-200
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	48kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human ZDHHC6
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Function: Palmitoylates calnexin (CALX), which is required for its association with the ribosome-translocon complex and efficient folding of glycosylated proteins.

Subcellular Location:

Membrane.

Similarity:

Belongs to the DHHC palmitoyltransferase family.
Contains 1 DHHC-type zinc finger.

Database links:

[Entrez Gene: 64429](#) Human

[Entrez Gene: 66980](#) Mouse

[Entrez Gene: 361771](#) Rat

[SwissProt: Q9H6R6](#) Human

[SwissProt: Q9CPV7](#) Mouse

[Unigene: 196990](#) Human

[Unigene: 386789](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.