

Rabbit Anti-MFAP3/AP Conjugated antibody

SL18823R-AP

Product Name	Anti-MFAP3/AP
Chinese Name	碱性磷酸酶（AP）标记的微丝相关蛋白 3 抗体
Alias	MFAP3; MFAP3_HUMAN; Microfibril-associated glycoprotein 3; Microfibrillar associated protein 3.
Research Area	Cell biology immunology Signal transduction Binding protein
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human(predicted:Mouse,Rat,Cow,Rabbit) WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	40kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human MFAP3
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	Function: Component of the elastin-associated microfibrils.
Product Detail	Subunit: Contains 1 Ig-like C2-type (immunoglobulin-like) domain. Subcellular Location: Cell membrane.

Post-translational modifications:

Glycosylated.

Similarity:

Contains 1 Ig-like C2-type (immunoglobulin-like) domain.

Database links:

[Entrez Gene: 4238](#) Human

[Omim: 600491](#) Human

[SwissProt: P55082](#) Human

[Unigene: 432818](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.