



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

## Rabbit Anti-ZNF576/Cy5.5 Conjugated antibody

SL19525R-Cy5. 5

<b>Product Name</b>	Anti-ZNF576/Cy5.5
<b>Chinese Name</b>	Cy5.5 标记的 Zinc finger protein576 抗体
<b>Alias</b>	FLJ22700; MGC2508; Zinc finger protein 576; ZN576_HUMAN; ZNF 576; ZNF576.
<b>Research Area</b>	transcriptional regulatory factor Zinc finger protein Epigenetics
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Human,Pig,Cow,Horse,Sheep) ICC/IF=1:50-200,IF=1:100-500
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	19kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human ZNF576
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Storage</b>	<b>Function:</b> May be involved in transcriptional regulation.
<b>Product Detail</b>	<b>Subcellular Location:</b> Nucleus.



**Similarity:**

Belongs to the krueppel C2H2-type zinc-finger protein family.  
Contains 4 C2H2-type zinc fingers.

**Database links:**

[Entrez Gene: 79177](#) Human

[SwissProt: Q9H609](#) Human

[Unigene: 11110](#) Human

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.