

Rabbit Anti-RASL12/AF350 Conjugated antibody

SL19738R-AF350

Product Name	Anti-RASL12/AF350
Chinese Name	AF350 标记的 RAS 样蛋白 12 抗体
Alias	Ras family member Ris; RAS like family 12 protein; RAS like protein 12; RAS like protein Ris; Ras-like protein family member 12; RAS-like protein Ris; RASL12; RASLC_HUMAN; RIS.
Research Area	Cell biology Signal transduction G protein signal
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Rat(predicted:Human,Mouse,Pig,Cow,Horse,Rabbit,Zebrafish) IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	30kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human RASL12
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Similarity: Belongs to the small GTPase superfamily. Ras family. Database links:



[Entrez Gene: 51285](#) Human

[Entrez Gene: 415530](#) Chicken

[Entrez Gene: 614252](#) Cow

[Entrez Gene: 70784](#) Mouse

[Entrez Gene: 315762](#) Rat

[SwissProt: Q08E00](#) Cow

[SwissProt: Q9NYN1](#) Human

[SwissProt: Q08AT1](#) Mouse

[Unigene: 27018](#) Human

[Unigene: 35827](#) Mouse

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.