

Rabbit Anti-SLC35D2/AP Conjugated antibody

SL19824R-AP

Product Name	Anti-SLC35D2/AP
Chinese Name	碱性磷酸酶（AP）标记的溶质载体家族蛋白 35 成员 D2 抗体
Alias	Fringe connection; HFRC 1; Hfrc; HFRC1; Homolog of Fringe connection protein 1; MGC117215; MGC142139; Nucleotide sugar transporter UGTrel8; S35D2_HUMAN; SLC35D2; Solute carrier family 35 member D2; SQV 7L; SQV7-like protein; SQV7L; UDP N acetylglucosamine transporter; UDP-galactose transporter-related protein 8; UDP-N-acetylglucosamine/UDP-glucose/GDP-mannose transporter; UGTrel 8; UGTrel8.
Research Area	Tumour Cell biology Signal transduction
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human
Applications	WB=1:500-2000 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	37kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human SLC35D2
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Product Detail	background:

Nucleotide sugars, which are synthesized in the cytosol or the nucleus, are high-energy donor substrates for glycosyltransferases located in the lumen of the endoplasmic reticulum and Golgi apparatus. Translocation of nucleotide sugars from the cytosol into the lumen compartment is mediated by specific nucleotide sugar transporters, such as SLC35D2 (Suda et al., 2004 [PubMed 15082721]).[supplied by OMIM, Mar 2008]

Function:

Antiporter transporting nucleotide sugars such as UDP-N-acetylglucosamine (UDP-GlcNAc), UDP-glucose (UDP-Glc) and GDP-mannose (GDP-Man) pooled in the cytosol into the lumen of the Golgi in exchange for the corresponding nucleosides monophosphates (UMP for UDP-sugars and GMP for GDP-sugars). May take part in heparan sulfate synthesis by supplying UDP-Glc-NAc, the donor substrate, and thus be involved in growth factor signaling.

Subcellular Location:

Golgi apparatus membrane.

Tissue Specificity:

Highly expressed in heart, kidney, small intestine, placenta, lung and peripheral blood leukocyte. Weakly expressed in skeletal muscle and spleen. Not expressed in brain, colon and thymus.

Similarity:

Belongs to the TPT transporter family. SLC35D subfamily.

Database links:

[Entrez Gene: 11046](#) Human

[Omim: 609182](#) Human

[SwissProt: Q76EJ3](#) Human

[Unigene: 494556](#) Human

[Unigene: 731624](#) Human

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com
