

## Rabbit Anti-PrV gE antibody

SL20521R

<b>Product Name</b>	PrV gE
<b>Chinese Name</b>	猪伪狂犬病毒 gE 抗体
<b>Alias</b>	Envelope glycoprotein E, gE, GE_SUHVR, Pseudorabies virus, Suid herpesvirus 1.
<b>Research Area</b>	Bacteria and viruses
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:PrV) ELISA=1:5000-10000
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Theoretical molecular weight</b>	62kDa
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human PrV gE: 441-540/577
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Buffer Solution</b>	(predicted:PrV)1M TBS(pH7.4) with 1% BSA, (predicted:PrV)3% Proclin300 and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
<b>Attention</b>	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
<b>PubMed</b>	<a href="#">PubMed</a>
<b>Product Detail</b>	