

Rabbit Anti-Glu-Glu Tag/AP Conjugated antibody

SL2115R-AP

Product Name	Anti-Glu-Glu Tag/AP
Chinese Name	碱性磷酸酶（AP）标记的 Glu-Glu tag 抗体
Alias	EYMPME epitope tag; EYMPME tag; Glu-Glu epitope tag.
Product Type	Tag anti
Research Area	Cell type markers
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Species independent WB=1:1000-5000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide of Glu-Glu Tag
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	background: Glu-Glu Tag is a commonly used epitope tag engineered onto the N- or C-terminus of a protein of interest so that the tagged protein can be analyzed and visualized.Because of the small size of the epitope, it is unlikely to alter the activity of the cloned sequence.

Important Note:



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.