

## Rabbit Anti-BAFFR/AF350 Conjugated antibody

SL2472R-AF350

<b>Product Name</b>	Anti-BAFFR/AF350
<b>Chinese Name</b>	AF350 标记的 B 细胞活化因子受体抗体
<b>Alias</b>	B Cell Activating Factor Receptor; B-cell-activating factor receptor; BAFF R; BAFF receptor; BAFF-R; BAFFR; BLyS receptor 3; BlySR3; BR 3; BR3; BROMIX; CD 268; CD268; CD268 antigen; CVID4; MGC138235; OTTHUMP00000028746; Prolixin; TNFRSF 13C; TNFRSF13C; TR13C_HUMAN; Tumor Necrosis Factor Receptor Subunit Member 13C; Tumor necrosis factor receptor superfamily member 13C.
<b>Research Area</b>	Cell biology Signal transduction Growth factors and hormones
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Human,Mouse,Rat,Dog,Cow,Horse,Rabbit) IF=1:100-500
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	19kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human BAFFR
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol
<b>Storage</b>	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Product Detail</b>	<b>background:</b> B cell-activating factor (BAFF) enhances B-cell survival in vitro and is a

regulator of the peripheral B-cell population. Overexpression of Baff in mice results in mature B-cell hyperplasia and symptoms of systemic lupus erythematosus (SLE). Also, some SLE patients have increased levels of BAFF in serum. Therefore, it has been proposed that abnormally high levels of BAFF may contribute to the pathogenesis of autoimmune diseases by enhancing the survival of autoreactive B cells. The protein encoded by this gene is a receptor for BAFF and is a type III transmembrane protein containing a single extracellular cysteine-rich domain. It is thought that this receptor is the principal receptor required for BAFF-mediated mature B-cell survival. [provided by RefSeq].

**Function:**

B-cell receptor specific for TNFSF13B/TALL1/BAFF/BLyS. Promotes the survival of mature B-cells and the B-cell response.

**Subcellular Location:**

Membrane; Single-pass type III membrane protein (Probable).

**Tissue Specificity:**

Highly expressed in spleen and lymph node, and in resting B-cells. Detected at lower levels in activated B-cells, resting CD4+ T-cells, in thymus and peripheral blood leukocytes.

**DISEASE:**

Defects in TNFRSF13C are the cause of immunodeficiency common variable type 4 (CVID4) [MIM:613494]; also called antibody deficiency due to BAFFR defect. CVID4 is a primary immunodeficiency characterized by antibody deficiency, hypogammaglobulinemia, recurrent bacterial infections and an inability to mount an antibody response to antigen. The defect results from a failure of B-cell differentiation and impaired secretion of immunoglobulins; the numbers of circulating B-cells is usually in the normal range, but can be low.

**Similarity:**

Contains 1 TNFR-Cys repeat.

**Database links:**

[Entrez Gene: 115650](#) Human

[Entrez Gene: 72049](#) Mouse

[Entrez Gene: 500910](#) Rat

[Omim: 606269](#) Human



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

---

[SwissProt: Q96RJ3](#) Human

[SwissProt: Q9D8D0](#) Mouse

[Unigene: 344088](#) Human

[Unigene: 240047](#) Mouse

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.