

Rabbit Anti-Sec14L3/Cy5.5 Conjugated antibody

SL4448R-Cy5. 5

Product Name	Anti-Sec14L3/Cy5.5
Chinese Name	Cy5.5 标记的 SEC14 样蛋白 L3 抗体
Alias	45 kDa secretory protein; MGC119366; MGC119367; MGC119368; S14L3_HUMAN; SEC14-like protein 3; Sec14I3; TAP2; Tocopherol-associated protein 2.
Research Area	Tumour Cell biology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human,Rat(predicted:Mouse) IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	46kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human Sec14L3
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	Function: Probable hydrophobic ligand-binding protein; may play a role in the transport of hydrophobic ligands like tocopherol, squalene and phospholipids. Similarity:



Contains 1 CRAL-TRIO domain.

Contains 1 GOLD domain.

Database links:

[Entrez Gene: 380683](#) Mouse

[Entrez Gene: 64543](#) Rat

[SwissProt: Q9Z1J8](#) Rat

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.