

Rabbit Anti-ABHD15/AP Conjugated antibody

SL7946R-AP

Product Name	Anti-ABHD15/AP
Chinese Name	碱性磷酸酶 (AP) 标记的水解酶结构域蛋白 15 抗体
Alias	ABH15_HUMAN; Abhd15; Abhydrolase domain containing protein 15; Abhydrolase domain-containing protein 15; hypothetical protein LOC116236; LOC116236; RGD1304598.
Research Area	Cell biology immunology
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Mouse,Rat(predicted:Human,Dog,Pig,Cow,Horse,Sheep) WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	49kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	KLH conjugated synthetic peptide derived from human ABHD15
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	Subcellular Location: Secreted (Potential).
Product Detail	Similarity: Belongs to the AB hydrolase superfamily. AB hydrolase 4 family. Database links:



SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

UniProtKB/Swiss-Prot: Q6UXT9.2

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.