

## Rabbit Anti-FAM46C/Cy5 Conjugated antibody

SL8198R-Cy5

<b>Product Name</b>	Anti-FAM46C/Cy5
<b>Chinese Name</b>	Cy5 标记的 FAM46C 蛋白抗体
<b>Alias</b>	Family with sequence similarity 46 member C; Protein FAM46C.
<b>Research Area</b>	
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	Rat(predicted:Human,Mouse,Chicken,Dog,Pig,Cow,Horse,Rabbit) IF=1:100-500
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	43kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human FAM46C
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Storage</b>	<b>background:</b> [FUNCTION] Seems to enhance replication of some viruses, including yellow fever virus, in response to type I interferon.
<b>Product Detail</b>	<b>Similarity:</b> Belongs to the FAM46 family.



**Database links:**

[Entrez Gene: 54855](#) Human

[Entrez Gene: 74645](#) Mouse

[Entrez Gene: 310721](#) Rat

[Olim: 613952](#) Human

[SwissProt: Q5VWP2](#) Human

[SwissProt: Q5SSF7](#) Mouse

[SwissProt: Q5XIV0](#) Rat

[Unigene: 356216](#) Human

[Unigene: 226900](#) Mouse

[Unigene: 162620](#) Rat

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.