

Rabbit Anti-DENND2D/Cy5 Conjugated antibody

SL9668R-Cy5

Product Name	Anti-DENND2D/Cy5
Chinese Name	Cy5 标记的 DENN 域内含蛋白 2D 抗体
Alias	DEN2D_HUMAN; DENN domain-containing protein 2D; DENN/MADD domain containing 2D; DENND2D.
Research Area	
Immunogen Species	Rabbit
Clonality	Polyclonal
React Species	Human(predicted:Mouse,Rat,Dog,Pig,Cow,Horse,Rabbit,Sheep) IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Molecular weight	54 kDa
Form	Lyophilized or Liquid
Concentration immunogen	1mg/ml KLH conjugated synthetic peptide derived from human DENND2D
Lsotype	IgG
Purification	affinity purified by Protein A
Storage Buffer	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Storage	
Product Detail	background: DENND2D is a 471 amino acid protein that contains a dDENN domain, a DENN domain, and a uDENN domain and exists as two isoforms as a result of alternative splicing. The DENND2D protein is thought to target to actin filaments and control Rab9-dependent trafficking of mannose-6-phosphate receptor to lysosomes. The gene encoding DENND2D maps to human

chromosome 1, the largest human chromosome which spans about 260 million base pairs and makes up 8% of the human genome. Other notable genes located on chromosome 1 include LMNA, which is associated with the rare aging disease Hutchinson-Gilford progeria, and the MUTYH gene, which is partially responsible for familial adenomatous polyposis. Stickler syndrome, Parkinsons, Gaucher disease and Usher syndrome.

Function:

Guanine nucleotide exchange factor (GEF) which may activate RAB9A and RAB9B. Promotes the exchange of GDP to GTP, converting inactive GDP-bound Rab proteins into their active GTP-bound form.

Similarity:

Contains 1 dDENN domain.

Contains 1 DENN domain.

Contains 1 uDENN domain.

Database links:

UniProtKB/Swiss-Prot: Q9H6A0.2

Important Note:

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.