

## Rabbit Anti-C1ORF111/Cy5 Conjugated antibody

SL9791R-Cy5

<b>Product Name</b>	Anti-C1ORF111/Cy5
<b>Chinese Name</b>	Cy5 标记的 1 号染色体开放阅读框 111 抗体
<b>Alias</b>	chromosome 1 open reading frame 111; HSD20; hypothetical protein LOC284680; CA111_HUMAN.
<b>Research Area</b>	Cell biology immunology
<b>Immunogen Species</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>React Species</b>	(predicted:Human,Mouse,Rat,Dog,Horse) IF=1:100-500
<b>Applications</b>	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
<b>Molecular weight</b>	29kDa
<b>Form</b>	Lyophilized or Liquid
<b>Concentration</b>	1mg/ml
<b>immunogen</b>	KLH conjugated synthetic peptide derived from human C1ORF111
<b>Lsotype</b>	IgG
<b>Purification</b>	affinity purified by Protein A
<b>Storage Buffer</b>	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 1M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
<b>Storage</b>	<b>background:</b> The function of the C1orf111 protein is currently unknown.
<b>Product Detail</b>	<b>Database links:</b> <a href="#">Entrez Gene: 284680</a> Human



SunLong Biotech Co.,LTD  
Tel: 0086-571-56623320 Fax:0086-571-56623318  
E-mail:sales@sunlongbiotech.com  
www.sunlongbiotech.com

---

[SwissProt: Q5TOL3](#) Human

**Important Note:**

This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.