



SUNLONG

SunLong Biotech Co.,LTD

Tel: 0086-571-56623320 Fax:0086-571-56623318

E-mail:sales@sunlongbiotech.com

www.sunlongbiotech.com

Mouse Anti-Digoxin antibody

SLM-33046M

| | |
|--------------------------|--|
| Product Name | Digoxin |
| Chinese Name | 地高辛单克隆抗体 |
| Alias | Lanoxin. |
| Product Type | Small molecule anti |
| Research Area | Cell biology immunology Drugs and Compounds |
| Immunogen Species | Mouse |
| Clonality | Monoclonal |
| Clone NO. | 6H7 |
| React Species | Species independent ELISA=1:1000-10000 |
| Applications | not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Form | Liquid |
| Concentration | 1mg/ml |
| immunogen | KLH conjugated Digoxin |
| Lsotype | IgG |
| Purification | affinity purified by Protein G |
| Buffer Solution | Species independent1M TBS(pH7.4) with 1% BSA, Species independent3% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Attention | This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. |
| PubMed | PubMed Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat). |
| Product Detail | |

SWISS:

N/A

CAS:
20830-75-5

Product
Picture

