

Mouse Anti-ASFV p30 antibody

SLM-41382M

Product Name	ASFV p30
Chinese Name	非洲猪瘟 p30 单克隆抗体
Alias	Phosphoprotein p30; Ba71V-93; P30_ASFB7; p30; Phosphoprotein p32; p32; CP204L
Immunogen Species	Mouse
Clonality	Monoclonal
React Species	(predicted:ASFV) ELISA=1:5000-10000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Theoretical molecular weight	23.5kDa
Form	Lyophilized or Liquid
Concentration	1mg/ml
immunogen	Recombined ASFV p30 protein: 1-204/204
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:ASFV)1M TBS(pH7.4) with 1% BSA, (predicted:ASFV)3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed Modifies the subcellular distribution of heterogeneous nuclear ribonucleoprotein K (HNRNPK) and may contribute to modulate HNRNPK functions related to processing and export of mRNAs during ASFV infection.
Product Detail	SWISS: P34204



SunLong Biotech Co.,LTD
Tel: 0086-571-56623320 Fax:0086-571-56623318
E-mail:sales@sunlongbiotech.com
www.sunlongbiotech.com

Gene ID:
22220322