



Mouse Anti-HBsAg antibody

SLM-41522M

Product Name	HBsAg
Chinese Name	乙肝病毒表面抗原单克隆抗体
Immunogen Species	Mouse
Clonality	Monoclonal
React Species	(predicted:HBV) WB=1:500-2000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Form	Liquid
immunogen	Recombinant Hepatitis B Virus
Lsotype	IgG
Purification	affinity purified by Protein A
Buffer Solution	(predicted:HBV)1M PBS(pH7.4)
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	