




Rabbit Anti-PCK2 antibody

SLM-60457R

Product Name	PCK2
Chinese Name	
Immunogen Species	Rabbit
Clonality	Monoclonal
Clone NO.	E6C11
React Species	Human,Mouse,Rat
Applications	WB=1:500-2000,Flow-Cyt=1:50-100,IHC-P=1:100-500,IHC-F=1:100-500,IF=1:100-500 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	The cell membrane
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG/Kappa
Purification	Affinity purified by Protein A
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	催化草酰乙酸（OAA）转化为磷酸烯醇式丙酮酸（PEP），这是代谢途径中的限速步骤，通过乳酸和柠檬酸循环产生的其他前体生成葡萄糖。

**Product
Picture**



Cell line: HeLa

Fixative: 100% Ice-cold methanol

Permeabilization: 0.1% TritonX-100

Primary Ab dilution: 1:50

Primary incubation condition: 4°C overnight

Secondary Ab: Goat Anti-Mouse IgG

Nuclear counter stain: DAPI (Blue)

Comment: Color green is the positive signal for SLM-60457R