

Mouse Anti-LPP antibody

SLM-60459M

Product Name	LPP
Chinese Name	
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	C12A9
React Species	(predicted: Human, Mouse, Rat,) WB=1:500-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	cytoplasmic The cell membrane
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG1/Kappa
Purification	Affinity purified by Protein G
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	可能在维持细胞形状和运动的细胞粘附部位发挥结构作用。除了这些结构功能外，它还可能与信号事件和基因转录激活有关。参与从细胞粘附位点到 The nucleus 的 Signal transduction，使可溶性因子和细胞-细胞粘附位点产生的信号成功整合。它也被认为是一种支架蛋白，不同的蛋白复合物在细胞质和 The nucleus 中组装在其上。