

Mouse Anti-NAPRT1 antibody

SLM-60483M

Product Name	NAPRT1
Chinese Name	
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	G6A3
React Species	(predicted: Human,) WB=1:500-1000
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	The cell membrane
Form	Liquid
Concentration	1mg/ml
Lsotype	Ig2A/Kappa
Purification	Affinity purified by Protein G
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	催化从烟酸生物合成 NAD 的第一步，即从烟酸和 5-磷酸-D-核糖 1-磷酸合成 β-烟酸-D-核糖核苷酸的 ATP 依赖性合成。通过其在 NAD 生物合成中的作用帮助防止细胞氧化应激。