



Mouse Anti-TFG antibody

SLM-60513M

| | |
|------------------------------|---|
| Product Name | TFG |
| Chinese Name | |
| Immunogen Species | Mouse |
| Clonality | Monoclonal |
| Clone NO. | E2F8 |
| React Species | (predicted: Human,) WB=1:1000-5000,ICC/IF=1:100-500 |
| Applications | not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user. |
| Cellular localization | The cell membrane |
| Form | Liquid |
| Concentration | 1mg/ml |
| Lsotype | IgG2B/Kappa |
| Purification | Affinity purified by Protein G |
| Buffer Solution | 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol. |
| Storage | Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. |
| Attention | This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications. |
| PubMed | PubMed |
| Product Detail | 在内质网及其相关微管的正常动态功能中起作用。从内质网到高尔基体的分泌性货物运输所必需的。 |