

## Mouse Anti-MRPL45 antibody

SLM-60517M

**Product Name** MRPL45

**Chinese Name**

**Immunogen Species** Mouse

**Clonality** Monoclonal

**Clone NO.** B5A2

**React Species** (predicted:Human)

WB=1:500-1000,IHC-P=1:100-500,IHC-F=1:100-500,ICC/IF=1:100-500,IF=1:100-500  
(Paraffin sections need antigen repair)

**Applications**

not yet tested in other applications.

optimal dilutions/concentrations should be determined by the end user.

**Cellular localization** The cell membrane

**Form** Liquid

**Concentration** 1mg/ml

**Lsotype** IgG1/Kappa

**Purification** Affinity purified by Protein G

**Buffer Solution** 1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.

**Storage** Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.

**Attention** This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.

**PubMed** [PubMed](#)

**Product Detail**

哺乳动物 Mitochondrion 核糖体蛋白质由核基因编码, 有助于 Mitochondrion 内的蛋白质合成。Mitochondrion 核糖体由一个小的 28S 亚单位和一个大的 39S 亚单位组成。MRPL45 是一种 39S 亚单位蛋白。