

Mouse Anti-SNAP antibody

SLM-60522M

Product Name	SNAP
Chinese Name	
Immunogen Species	Mouse
Clonality	Monoclonal
Clone NO.	B9C1
React Species	(predicted: Human, Mouse, Rat,) WB=1:500-1000,ICC/IF=1:100-500
Applications	not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Cellular localization	cytoplasmic
Form	Liquid
Concentration	1mg/ml
Lsotype	IgG2B/lambda
Purification	Affinity purified by Protein G
Buffer Solution	1M TBS(pH7.4) with 1% BSA, 3% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Attention	This product as supplied is intended for research use only, not for use in human, therapeutic or diagnostic applications.
PubMed	PubMed
Product Detail	内质网和高尔基体之间的囊泡运输所需。与 GNA12 一起促进 CDH5 在质膜上的定位。